





Particular Sciences Ltd

Spinvette:— single step centrifugal DLS sample cleanup

Material characterisation technology



Dynamic Light Scattering measurement of real world samples often face the issue of tiny quantities of aggregated or large sized contaminant making the result difficult to interpret.



Coarse aggregate can be filtered with expensive nanofilters or centrifuged and spun to the bottom of the sample, then decanted or pipetted. Or you can use a Spinvette for single step clean-up and measure direct in the Zetasizer Nano. Save time and wasted sample.

A Spinvette adaptor to fit your Zetasizer Nano S or ZS or S90- complete with a starter set of 10 spinvettes is on offer for the special price of €725.00 vat extra (includes postal delivery).

Particular Sciences Ltd

2 Birch House
Rosemount Business Park
Ballycoolin Road
Dublin 11, Ireland
Phone: (353) 01 8205395
Fax: (353) 01 8228813
E-mail: info@particular.ie

- ⇒ Spinvettes are optical quality plastic and re-useable.
- ⇒ Work at 90 degrees, NIBS angles or forward scattering angles.
- ⇒ Small volumes easy to test.
- ⇒ Individual Spinvette adaptors to fit other makes of DLS available for €799.00 vat extra.

Spinvette is a product of Cusop Technology Inc.